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AGRICULTURE AND AGRIFOOD DEPARTMENT

MINISTRY OF PRIMARY RESOURCES AND TOURISM | BRUNEI DARUSSALAM

GUIDELINES FOR SAMPLE COLLECTION AND SUBMISSION

SALMONELLA REFERENCE LABORATORY (SAL)

Veterinary Laboratory Services Unit
Division of Livestock Industry
Department of Agriculture and Agrifood

First Edition 2020

Compiled by

Veterinary Laboratory Services Unit

Division of Livestock Industry

Department of Agriculture and Agrifood

Ministry of Primary Resources and Tourism

Brunei Darussalam

Copyright 2020

All rights reserved. No part of this book may be reproduced in any form or by any means without prior permission from the Veterinary Laboratory Services Unit, Division of Livestock Industry, Department of Agriculture and Agrifood, Ministry of Primary Resources and Tourism, Brunei Darussalam.

Perpustakaan Dewan Bahasa dan Pustaka Brunei Cataloguing-in-publication

BRUNEI Darussalam, Veterinary Laboratory Services Unit

Guidelines for Sample Collection and Submission – Salmonella Reference Laboratory /

Brunei Darussalam, Veterinary Laboratory Services Unit - Bandar Seri Begawan: Jabatan Pertanian dan Agrimakanan, 2020.

p. cm.

JP/BPPA/01/2020 (Kulit Keras)

1. Salmonella Reference I. Title.

636.089 BRU (DDC22)

TABLE OF CONTENTS

FOREWORD	1
INTRODUCTION	2
SAMPLES FOR SALMONELLA REFERENCE LABORATORY ANALYSIS	2
SAMPLE RETENTION TIME	4
CLIENT CHARTER	4
NEW CLIENTS	4

FOREWORD

Alhamdulillah, this Guideline for Sample Collection and Submission was developed to provide farmers, importers/ exporters of livestock and veterinary products, with information on sample collection and submission, type of samples and quantity of samples required for specific testing at Salmonella Reference Laboratory in Veterinary Laboratory Services Unit (VLS), Division of Livestock Industry, under the Department of Agriculture and Agrifood.

The detailed procedures and requirements in this booklet are based on the ISO/IEC 17025 and test method procedures for serotyping and antimicrobial susceptibility testing requirement. This booklet also contains the client charters of each testing provided by VLS for client reference.

We hope that the Guideline for Sample Collection and Submission Salmonella Reference Laboratory (SAL) booklet can be used as a guide for farmers, importers/ exporters and the public who wish to send sample for serotyping and antimicrobial susceptibility testing in VLS.

Finally, I would like to applaud and congratulate the VLS unit for their great effort in publishing this guideline. May Allah Subhanahu wa Ta'ala bless you all.

Wassalam

A handwritten signature in blue ink, consisting of several overlapping loops and a long horizontal stroke extending to the right.

[HAJAH MASLIANA YULIAH BINTI ABDULLAH]

Acting Director of Agriculture and Agrifood

INTRODUCTION

Salmonella Reference Laboratory under Veterinary Laboratory Services (VLS) offers the following analysis:

- (i) Salmonella serotyping
- (ii) Antimicrobial Susceptibility Testing

More than 2,500 genus *Salmonella* are reported worldwide and validated by World Health of Organisation (WHO). Serotyping of *salmonella* species is aimed for traceability of diseases especially in the farm so that outbreak could be prevented and controlled in an earlier manner. Knowing that certain *salmonella* serotypes can cause massive production losses.

The genus *salmonella* is then tested for antimicrobial susceptibility. The objective is to identify the prevalent use of antibiotics in the farm. This surveillance activity is to protect the efficiency of drugs and delay development of its drug resistance.

The information in this guideline is prepared to enable you to understand, develop a sampling plan for your operation prior to testing, sampling and sending the samples to VLS.

SAMPLES FOR SALMONELLA REFERENCE LABORATORY ANALYSIS

The samples tested for serological testing in Salmonella Reference Laboratory is in a form of pure bacterial culture. These are initially depicted from isolated *salmonella* from diagnostics samples as well as veterinary products samples submitted to Veterinary Laboratory Services through Veterinary Microbiology Laboratory and Microbiology Laboratory.

The samples such as **fresh faeces, drag swabs, cloacal swabs, boot swabs, litter, water, animal feed, infertile eggs, DOC box liners, organs, meat (chicken, duck, ostrich, cattle, buffalo, lamb, etc.), processed products and table eggs** are tested first either through Veterinary Microbiology Laboratory of Animal Diagnostic Sub-unit or

Microbiology Laboratory of Veterinary Public Health Sub-unit and continue in the Salmonella Reference Laboratory for confirmation of bacterial culture serotyping. Overall, these are aimed to support Department's National Salmonella Program and surveillance against antimicrobial resistance.

Clients are advised to observe the samples requirements for different sample types and delivery time stipulated in **Guidelines for Sample Collection and Submission Animal Diagnostic (ADG)** and **Guidelines for Sample Collection and Submission Veterinary Public Health (VPH)**.

The following procedures are usually employed in Salmonella Reference Laboratory:

- a) Upon receiving the pure bacterial culture samples, all existing clients have agreed that any positive samples tested are required for further *salmonella* serological testing as well as antimicrobial susceptibility testing, where relevant.
- b) As soon as the laboratory isolated the positive *salmonella*, Veterinary Microbiology Laboratory (Animal Diagnostic Sub-unit) or Microbiology Laboratory (Veterinary Public Health Sub-unit) technician will, on behalf of the clients, fill in a Salmonella Reference Test Request Form, VLS-SOP-09-06 (stating the name of the client) which is readily provided at the counter. Make sure the form is completed with all appropriate and necessary information required.
- c) For information, the Veterinary Microbiology Laboratory (Animal Diagnostic Sub-unit) or Microbiology Laboratory (Veterinary Public Health Sub-unit) technician will transport the pure bacterial culture in a highly sealed manner, keep in biohazard bag and in a closed container.
- d) The culture is required to be transported immediately to Salmonella Reference Laboratory within 12-24 hours.
- e) All accompanying **Salmonella Reference Test Request Form** shall be completely filled with **BLUE** ink **ONLY**.

SAMPLE RETENTION TIME

All samples received will be used for the whole process of testing. Salmonella Reference Laboratory may keep the growth culture frozen at -80°C by preserving them using Microbank or using glycerol-saline solution for collection culture and reference. Any cultures found impure or contaminated will be discarded immediately.

CLIENT CHARTER

The table below indicates the turnaround time required for each tests.

Table 1: Client charter

No.	Analysis activities/ Scope of Analysis	No. of working days
1	To Perform <i>Salmonella</i> Serotyping Analysis	5
2	To Perform Antimicrobial Susceptibility Testing	10

VLS always aim to produce reports as soon as it is ready.

The clients are encouraged to be aware of the turnaround time of the analysis for collecting test report by referring to the Case Reference Number. The test report must be collected by the clients themselves at the Veterinary Laboratory Services Unit reception counter.

NEW CLIENTS

New customer is required to fill up the New Testing Request Form, VLS-SOP-09-02 and attend a short briefing on the tests offered by VLS. Existing customers requesting for new test are also required to fill the above form. The form is available upon request at the reception counter.

**FOR FURTHER ENQUIRIES
PLEASE CALL / CONTACT
VETERINARY LABORATORY SERVICES UNIT**

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